


How to use tickatlab_anishare_libreoffice_macro with LibreOffice

Hint: Tested with LibreOffice Version 6.0.4.2.


Please contact Fabian Monheim (fabian.monheim@leibniz-fli.de , 03641-65-6872) if you have a question or the macro is not working properly.

1. Download zip with macro and unzip it
2. Open LibreOffice Calc
3. Go to *Tools* → *Macros* → *Organize Macros* → *LibreOffice Basic ...*




4. Click the *Organizer...* button on the right side

5. Go to *Libraries*



6. *Import...* the library with the macro. Choose at the import dialog the *script.xlb* file.


7. Click *OK*



Now the Marco TickAtLab_AnShare is imported and can be executed e.g. on this way: Tools → Macros → Run Macro... → My Macros → TickAtLab → TickAtLab → TickAtLab_AnShare
To simplify the execution it is possible to create button at the toolbar to run this macro and only one click is necessary to execute the macro.

Create Button at the toolbar to run the macro:

1. Go to *Tools* → *Customize...*



2. Choose *Toolbars*

3. Click on *All commands*


3.1. Choose *Macros* (second last choice)

4. Select Function: *My Macros* → *TickAtLab* → *TickAtLab* → *TickAtLab_AnShare*

4.1. Click the right arrow button

4.2. Modify the place of the *TickAtLab_AnShare* button at the toolbar with the up and down arrows on the right side.

4.3. Click the *OK* button



5. *TickAtLab_AniShare* button is ready to use

